

Product datasheet for **GA112340**

Human TUG (ASPSCR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ASPSCR1 CRISPRa kit - CRISPR gene activation of human ASPSCR1 tether for SLC2A4, UBX domain containing
Symbol:	ASPSCR1
UniProt ID:	Q9BZE9
Synonyms:	ASPCR1; ASPL; ASPS; RCC17; TUG; UBXD9; UBXN9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)