


Product datasheet for **GA112341**

Human TMEM223 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM223 CRISPRa kit - CRISPR gene activation of human transmembrane protein 223
Symbol:	TMEM223
Components:	GA112341G1, TMEM223 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112341G2, TMEM223 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112341G3, TMEM223 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A0PIW6
Synonyms:	MGC3196
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)