

Product datasheet for **GA112347**

Human FASTKD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FASTKD3 CRISPRa kit - CRISPR gene activation of human FAST kinase domains 3
Symbol:	FASTKD3
Components:	GA112347G1, FASTKD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112347G2, FASTKD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112347G3, FASTKD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14CZ7</u>
Synonyms:	FLJ23274; MGC5297; MGC142123
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)