

Product datasheet for **GA112348**

Human TMEM109 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM109 CRISPRa kit - CRISPR gene activation of human transmembrane protein 109
Symbol:	TMEM109
Components:	GA112348G1, TMEM109 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112348G2, TMEM109 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112348G3, TMEM109 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BVC6</u>
Synonyms:	MGC5508
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)