

Product datasheet for **GA112369**

Human SARG (C1orf116) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C1orf116 CRISPRa kit - CRISPR gene activation of human chromosome 1 open reading frame 116
Symbol:	C1orf116
Components:	GA112369G1, SARG gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112369G2, SARG gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112369G3, SARG gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BW04</u>
Synonyms:	SARG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)