

Product datasheet for **GA112370**

Human TAFID activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAFID CRISPRa kit – CRISPR gene activation of human TATA-box binding protein associated factor, RNA polymerase I subunit D
Symbol:	TAFID
Components:	GA112370G1, TAFID gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112370G2, TAFID gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112370G3, TAFID gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H5J8</u>
Synonyms:	JOSD3; RAFI41; TAF(I)41; TAFI41
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

