

Product datasheet for **GA112380**

Human PHF23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHF23 CRISPRa kit - CRISPR gene activation of human PHD finger protein 23
Symbol:	PHF23
UniProt ID:	<u>Q9BUL5</u>
Synonyms:	hJUNE-1b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)