


Product datasheet for **GA112388**

Human GDPD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GDPD3 CRISPRa kit - CRISPR gene activation of human glycerophosphodiester phosphodiesterase domain containing 3
Symbol:	GDPD3
Components:	GA112388G1, GDPD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112388G2, GDPD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112388G3, GDPD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7L5L3
Synonyms:	FLJ22603; glycerophosphodiester phosphodiesterase domain containing 3; MGC4171; MGC4171, FLJ22603
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)