

## Product datasheet for **GA112395**

### Human **C7orf23 (TMEM243) activation kit by CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM243 CRISPRa kit - CRISPR gene activation of human transmembrane protein 243
Symbol:	TMEM243
UniProt ID:	<a href="#">Q9BU79</a>
Synonyms:	C7orf23; MM-TRAG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)