

Product datasheet for **GA112397**

Human LILRA6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LILRA6 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor A6
Symbol:	LILRA6
Components:	GA112397G1, LILRA6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112397G2, LILRA6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112397G3, LILRA6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	CD85b; ILT-8; ILT5; ILT8; LILRB3; LILRB6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)