

## Product datasheet for **GA112404**

### Human ZNF343 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ZNF343 CRISPRa kit - CRISPR gene activation of human zinc finger protein 343
Symbol:	ZNF343
UniProt ID:	<a href="#">Q6P1L6</a>
Synonyms:	dj734P14.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)