

Product datasheet for **GA112412**

Human TMEM43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM43 CRISPRa kit - CRISPR gene activation of human transmembrane protein 43
Symbol:	TMEM43
UniProt ID:	Q9BTV4
Synonyms:	ARVC5; ARVD5; EDMD7; LUMA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)