

Product datasheet for **GA112413**

Human IRX6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IRX6 CRISPRa kit - CRISPR gene activation of human iroquois homeobox 6
Symbol:	IRX6
Components:	GA112413G1, IRX6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112413G2, IRX6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112413G3, IRX6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P78412
Synonyms:	IRX-3; IRX7; IRXB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)