

Product datasheet for **GA112414**

Human Iroquois homeobox protein 3 (IRX3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IRX3 CRISPRa kit - CRISPR gene activation of human iroquois homeobox 3
Symbol:	IRX3
UniProt ID:	P78415
Synonyms:	IRX-1; IRXB1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)