

Product datasheet for **GA112441**

Human LRFN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRFN3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat and fibronectin type III domain containing 3
Symbol:	LRFN3
UniProt ID:	Q9BTN0
Synonyms:	FIGLER1; SALM4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)