

## Product datasheet for **GA112443**

### Human HAUS3 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                            |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| <b>Description:</b>  | HAUS3 CRISPRa kit - CRISPR gene activation of human HAUS augmin like complex subunit 3                                                                                                                                                                                                                     |
| <b>Symbol:</b>       | HAUS3                                                                                                                                                                                                                                                                                                      |
| <b>Components:</b>   | GA112443G1, HAUS3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112443G2, HAUS3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112443G3, HAUS3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q68CZ6</a></u>                                                                                                                                                                                                                                                                              |
| <b>Synonyms:</b>     | C4orf15; dgt3; IT1                                                                                                                                                                                                                                                                                         |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:

