


Product datasheet for **GA112447**

Human WDR25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR25 CRISPRa kit - CRISPR gene activation of human WD repeat domain 25
Symbol:	WDR25
Components:	GA112447G1, WDR25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112447G2, WDR25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112447G3, WDR25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q64LD2</u>
Synonyms:	C14orf67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)