

Product datasheet for **GA112448**

Human PA1 (PAGR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAGR1 CRISPRa kit - CRISPR gene activation of human PAXIP1 associated glutamate rich protein 1
Symbol:	PAGR1
Components:	GA112448G1, PA1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112448G2, PA1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112448G3, PA1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BTK6</u>
Synonyms:	C16orf53; GAS; PA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)