

Product datasheet for **GA112464**

Human EPS8 like protein 3 (EPS8L3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EPS8L3 CRISPRa kit - CRISPR gene activation of human EPS8 like 3
Symbol:	EPS8L3
UniProt ID:	Q8TE67
Synonyms:	EPS8R3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)