

Product datasheet for **GA112469**

Human SPAG16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPAG16 CRISPRa kit - CRISPR gene activation of human sperm associated antigen 16
Symbol:	SPAG16
UniProt ID:	Q8N0X2
Synonyms:	PF20; WDR29
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)