

Product datasheet for **GA112470**

Human TMEM231 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM231 CRISPRa kit - CRISPR gene activation of human transmembrane protein 231
Symbol:	TMEM231
Components:	GA112470G1, TMEM231 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112470G2, TMEM231 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112470G3, TMEM231 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6L2</u>
Synonyms:	ALYE870; JBTS20; MKS11; PRO1886
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)