

Product datasheet for **GA112476**

Human ARMH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARMH3 CRISPRa kit – CRISPR gene activation of human armadillo like helical domain containing 3
Symbol:	ARMH3
Components:	GA112476G1, ARMH3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112476G2, ARMH3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112476G3, ARMH3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T2E6</u>
Synonyms:	RP11-190J1.9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

