

Product datasheet for **GA112477**

Human E3 ubiquitin protein ligase MUL1 (MUL1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MUL1 CRISPRa kit - CRISPR gene activation of human mitochondrial E3 ubiquitin protein ligase 1
Symbol:	MUL1
UniProt ID:	Q969V5
Synonyms:	C1orf166; GIDE; MAPL; MULAN; RNF218
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)