

Product datasheet for **GA112486**

Human **RIC3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RIC3 CRISPRa kit - CRISPR gene activation of human RIC3 acetylcholine receptor chaperone
Symbol:	RIC3
UniProt ID:	<u>Q7Z5B4</u>
Synonyms:	AYST720; PRO1385
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)