


## Product datasheet for **GA112494**

### Human SNRNP25 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SNRNP25 CRISPRa kit - CRISPR gene activation of human small nuclear ribonucleoprotein U11/U12 subunit 25
Symbol:	SNRNP25
Components:	GA112494G1, SNRNP25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112494G2, SNRNP25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112494G3, SNRNP25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q9BV90</a></u>
Synonyms:	C16orf33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)