

Product datasheet for **GA112500**

Human SH3TC2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3TC2 CRISPRa kit - CRISPR gene activation of human SH3 domain and tetratricopeptide repeats 2
Symbol:	SH3TC2
UniProt ID:	Q8TF17
Synonyms:	CMT4C; MNMN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)