

Product datasheet for **GA112516**

Human AKIRIN1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | AKIRIN1 CRISPRa kit - CRISPR gene activation of human akirin 1 |
| Symbol: | AKIRIN1 |
| UniProt ID: | Q9H9L7 |
| Synonyms: | C1orf108; STRF2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)