

Product datasheet for **GA112531**

Human DHX40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DHX40 CRISPRa kit - CRISPR gene activation of human DEAH-box helicase 40
Symbol:	DHX40
Components:	GA112531G1, DHX40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112531G2, DHX40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112531G3, DHX40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IX18</u>
Synonyms:	ARG147; DDX40; PAD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)