

Product datasheet for **GA112538**

Human ZNF329 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF329 CRISPRa kit - CRISPR gene activation of human zinc finger protein 329
Symbol:	ZNF329
Components:	GA112538G1, ZNF329 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112538G2, ZNF329 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112538G3, ZNF329 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UD4</u>
Synonyms:	FLJ12586
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)