

Product datasheet for **GA112551**

Human QTRTD1 (QTRT2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	QTRT2 CRISPRa kit - CRISPR gene activation of human queuine tRNA-ribosyltransferase accessory subunit 2
Symbol:	QTRT2
Components:	GA112551G1, QTRTD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112551G2, QTRTD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112551G3, QTRTD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H974</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)