


Product datasheet for **GA112553**

Human YRDC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	YRDC CRISPRa kit - CRISPR gene activation of human yrdC N6-threonylcarbamoyltransferase domain containing
Symbol:	YRDC
Components:	GA112553G1, YRDC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112553G2, YRDC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112553G3, YRDC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86U90</u>
Synonyms:	DRIP3; IRIP; SUA5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)