

Product datasheet for **GA112563**

Human PRKRIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKRIP1 CRISPRa kit - CRISPR gene activation of human PRKR interacting protein 1
Symbol:	PRKRIP1
Components:	GA112563G1, PRKRIP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112563G2, PRKRIP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112563G3, PRKRIP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H875
Synonyms:	C114; KRBOX3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)