


Product datasheet for **GA112580**

Human WDR59 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR59 CRISPRa kit - CRISPR gene activation of human WD repeat domain 59
Symbol:	WDR59
Components:	GA112580G1, WDR59 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112580G2, WDR59 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112580G3, WDR59 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6PJI9
Synonyms:	FP977
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)