

Product datasheet for **GA112599**

Human LMAN1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LMAN1L CRISPRa kit - CRISPR gene activation of human lectin, mannose binding 1 like
Symbol:	LMAN1L
Components:	GA112599G1, LMAN1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112599G2, LMAN1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112599G3, LMAN1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HAT1</u>
Synonyms:	ERGIC-53L; ERGL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)