

Product datasheet for **GA112600**

Human ZNF385D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF385D CRISPRa kit - CRISPR gene activation of human zinc finger protein 385D
Symbol:	ZNF385D
Components:	GA112600G1, ZNF385D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112600G2, ZNF385D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112600G3, ZNF385D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6B1</u>
Synonyms:	ZNF659
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)