

Product datasheet for **GA112649**

Human **MOBKL2B (MOB3B)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MOB3B CRISPRa kit - CRISPR gene activation of human MOB kinase activator 3B
Symbol:	MOB3B
UniProt ID:	Q86TA1
Synonyms:	C9orf35; MOB1D; MOBKL2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)