

Product datasheet for **GA112664**

Human LONRF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LONRF3 CRISPRa kit - CRISPR gene activation of human LON peptidase N-terminal domain and ring finger 3
Symbol:	LONRF3
UniProt ID:	Q496Y0
Synonyms:	RNF127
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)