

## Product datasheet for **GA112665**

### Human PIP4K2C activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIP4K2C CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-5-phosphate 4-kinase type 2 gamma
Symbol:	PIP4K2C
UniProt ID:	<a href="#">Q8TBX8</a>
Synonyms:	PIP5K2C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

