

Product datasheet for **GA112666**

Human TMC5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMC5 CRISPRa kit - CRISPR gene activation of human transmembrane channel like 5
Symbol:	TMC5
Components:	GA112666G1, TMC5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112666G2, TMC5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112666G3, TMC5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UXY8</u>
Synonyms:	FLJ13593
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)