

Product datasheet for **GA112677**

Human PDZD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDZD3 CRISPRa kit - CRISPR gene activation of human PDZ domain containing 3
Symbol:	PDZD3
Components:	GA112677G1, PDZD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112677G2, PDZD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112677G3, PDZD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UT5</u>
Synonyms:	IKEPP; NHERF4; PDZK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)