

Product datasheet for **GA112680**

Human **TM4SF20** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TM4SF20 CRISPRa kit - CRISPR gene activation of human transmembrane 4 L six family member 20
Symbol:	TM4SF20
UniProt ID:	Q53R12
Synonyms:	PRO994; SLI5; TCCE518
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

