


Product datasheet for **GA112697**

Human THSD4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	THSD4 CRISPRa kit - CRISPR gene activation of human thrombospondin type 1 domain containing 4
Symbol:	THSD4
Components:	GA112697G1, THSD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112697G2, THSD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112697G3, THSD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZMP0</u>
Synonyms:	ADAMTSL-6; ADAMTSL6; FVSY9334; PRO34005
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)