

## Product datasheet for **GA112737**

### Human DOK3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DOK3 CRISPRa kit - CRISPR gene activation of human docking protein 3
Symbol:	DOK3
UniProt ID:	<a href="#">Q7L591</a>
Synonyms:	DOKL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)