

## Product datasheet for **GA112738**

### Human ABIN3 (TNIP3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNIP3 CRISPRa kit - CRISPR gene activation of human TNFAIP3 interacting protein 3
Symbol:	TNIP3
UniProt ID:	<a href="#">Q96KP6</a>
Synonyms:	ABIN-3; LIND
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)