

Product datasheet for **GA112739**

Human KIAA0319L activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | KIAA0319L CRISPRa kit - CRISPR gene activation of human KIAA0319 like |
| Symbol: | KIAA0319L |
| UniProt ID: | <u>Q8IZA0</u> |
| Synonyms: | FLJ44532; KIAA1837 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)