

Product datasheet for **GA112744**

Human SLC35E1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35E1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member E1
Symbol:	SLC35E1
Components:	GA112744G1, SLC35E1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112744G2, SLC35E1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112744G3, SLC35E1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96K37
Synonyms:	DKFZp564G0462; FLJ14251; FLJ36689; MGC44954
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)