

Product datasheet for **GA112749**

Human LPPR3 (PLPPR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPPR3 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase related 3
Symbol:	PLPPR3
UniProt ID:	Q6T4P5
Synonyms:	LPPR3; LPR3; PRG-2; PRG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)