

## Product datasheet for **GA112750**

### Human C10orf81 (PLEKHS1) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | PLEKHS1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing S1                                                                                                                                                                                                                      |
| Symbol:       | PLEKHS1                                                                                                                                                                                                                                                                                                             |
| Components:   | GA112750G1, C10orf81 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112750G2, C10orf81 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112750G3, C10orf81 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q5SXH7</a>                                                                                                                                                                                                                                                                                              |
| Synonyms:     | bA211N11.2; C10orf81; HEL185                                                                                                                                                                                                                                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)