

Product datasheet for **GA112775**

Human PLEKHH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHH3 CRISPRa kit - CRISPR gene activation of human pleckstrin homology, MyTH4 and FERM domain containing H3
Symbol:	PLEKHH3
Components:	GA112775G1, PLEKHH3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112775G2, PLEKHH3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112775G3, PLEKHH3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z736</u>
Synonyms:	FLJ21019
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)