

Product datasheet for **GA112794**

Human FOXRED2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FOXRED2 CRISPRa kit - CRISPR gene activation of human FAD dependent oxidoreductase domain containing 2
Symbol:	FOXRED2
Components:	GA112794G1, FOXRED2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112794G2, FOXRED2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112794G3, FOXRED2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IWF2</u>
Synonyms:	ERFAD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)